



CE  
2797  
EC REP  
Emergo Europe  
Prinsessegracht 20  
2514 AP  
The Hague  
The Netherlands



QuikClot EMS Rolled Gauze

REF 475

**QuikClot EMS Rolled Gauze** is a sterile, soft white 3 in. x 4 ft. nonwoven hemostatic gauze impregnated with kaolin, an inorganic mineral that activates Factor XII,<sup>1</sup> which in turn accelerates the body's natural clotting ability.<sup>2,3,4,5</sup> Allowing you to control bleeding fast,<sup>9</sup> potentially giving you more time to attend to other tasks as you move patients through the continuum of care. Compared to standard rolled gauze or other available hemostatic devices, QuikClot has superior results that indicate it can reduce blood loss,<sup>9</sup> and promote robust clot formation.

**SCAN TO PAGE**



**ARIS IS A DISTRIBUTOR FOR TELEFLEX PRODUCTS**

Write to us for a Quote: [productdev@aris-im.com](mailto:productdev@aris-im.com)

Contact: 96550253 / [6283 0939](tel:62830939)

Address: 18 howard Road #08-02, Singapore 369585