



# FOXSEAL™ CHEST SEAL


(EXCLUSIVE DISTRIBUTOR)

FOXSEAL™ Chest Seal is an innovative occlusive dressing designed for managing open chest wounds to help prevent pneumothorax.



- Dimensions: 13 x 13.3cm chest seal dressings
- Provided in one sterile foil pouch measuring 15.8 x 15.56in
- Valve-free.
- Not made with Latex.
- Durable packaging has 4 tear notches for opening.

**Provides dependable adhesion in the most challenging situations.**<sup>22</sup>



<<< More Product Info  
CONTACT US:  
productdev@aris-im.com or +65 96550253