





# NASOPHARYNGEAL AIRWAY



Part No	ID (mm)	OD (mm)	Length (mm)	Box
NS-12FR-GA-DT	1.5	4.0	89	10
NS-14FR-GA-DT	2.0	4.7	89	10
NS-16FR-GA-DT	2.5	5.6	89	10
NS-18FR-GA-DT	3.2	6.0	89	10
NS-20FR-GA-DT	3.9	6.7	114	10
NS-22FR-GA-DT	4.8	7.3	146	10
NS-24FR-GA-DT	5.5	8.0	159	10
NS-26FR-GA-DT	6.3	8.7	159	10
NS-28FR-GA-DT	6.8	9.3	159	10
NS-30FR-GA-DT	7.0	10.00	159	10
NS-32FR-GA-DT	7.6	10.7	159	10
NS-34FR-GA-DT	8.6	11.3	159	10
NS-36FR-GA-DT	9.0	12.00	159	10

The Robertazzi Nasopharyngeal Airway is designed to be soft and flexible to minimize discomfort to the patient, made from high quality neoprene material. It comes with smoothly rounded edges to provide additional comfort to the patient.

