



R-LIFT by MATBOCK



Product Specs:

Material: Ghost Heavy

Weights:

Weight with Combat Carbon Poles & Pack:
6.5 lbs (2.94 kgs)

Weight with Combat Carbon Poles:
6 lbs (2.72 kgs)

Dimensions:

- 90.5" L x 27" W

(229.87cm L x 68.58cm W) - Deployed

- 21.25" L x 15" W x 5" D

(53.97cm L x 38.1cm W x 12.7cm D) - Packed



SCAN ME