



# Nasopharyngeal Airway

- Features:
- Self lubricating
  - Soft
  - Economical
  - Adjustable flange
  - Sterile
  - Disposable



Part No	ID (mm)	OD (mm)	Length (mm)	Box
NS-12FR-GA-DT	1.5	4.0	89	10
NS-14FR-GA-DT	2.0	4.7	89	10
NS-16FR-GA-DT	2.5	5.6	89	10
NS-18FR-GA-DT	3.2	6.0	89	10
NS-20FR-GA-DT	3.9	6.7	114	10
NS-22FR-GA-DT	4.8	7.3	146	10
NS-24FR-GA-DT	5.5	8.0	159	10
NS-26FR-GA-DT	6.3	8.7	159	10
NS-28FR-GA-DT	6.8	9.3	159	10
NS-30FR-GA-DT	7.0	10.00	159	10
NS-32FR-GA-DT	7.6	10.7	159	10
NS-34FR-GA-DT	8.6	11.3	159	10
NS-36FR-GA-DT	9.0	12.00	159	10

SCAN TO PAGE



Write to us for a Quote: [productdev@aris-im.com](mailto:productdev@aris-im.com)  
 Contact: 96550253 / [6283 0939](tel:62830939)  
 Address: 18 Howard Road #08-02, Singapore 369585