




IV Splints

I.V.Board is made of thin aluminum plate padded with soft foam. It's flexible and can be used for stabilizing the wrist during vascular access procedures to prevent patients from flexing their wrist.

- Features:
1. Malleable for wrist part of patients
 2. Soft with foam cover
 3. No skin irritating and latex free
 4. Well extending and stabilizing the wrist
 5. Reusable and can be reshaped
 6. Comfortable and safe
 7. Easy to apply and remove
 8. Waterproof type
 9. Normal color is blue & white, others are

EMAIL TO US FOR DIMENSIONS:



Item Number			Length(mm)	Width(mm)	Thickness (mm)
AZ-ABF-01	Small	Neonate	67	26(18)	8
AZ-ABF-02	Medium	Infant	105	35(30)	8
AZ-ABF-03	Large	Pediatrics	138	50(37)	8
AZ-ABF-04	Extra Large	Adolescent	175	55	8
AZ-ABF-05	Extra Extra Large	Adult	260	90	8

SCAN TO PAGE



Write to us for a Quote: productdev@aris-im.com

Contact: 96550253 / 6283 0939

Address: 18 howard Road #08-02, Singapore 369585